

## **Covid Product profile**

 Product name + short description (advantage of product in treatment or prevention of Covid)
 3-5 sentences:

Products against encapsulated virus/Covid 19

Personal care products:

INNO-SEPT cleaning and disinfection in one step
LILLA-SEPT cleaning and disinfection in one step
INNO-SEPT GEL EXTRA alcohol based hand disinfection
INNO-SEPT FRESH DISINFECTANT SOLUTION hand and surface
ready-to-use disinfection solution
INNO-SEPT FRESH DISINFECTANT WIPES hand and surface
disinfection wipes

2. link/webpage of product or company:

www.innoveng1.hu

3. contact:

Anita Tóth

anita.toth@innoveng1.hu

Kamran Afrid/CEO

kamran.afrid@innoveng1.hu

- 4. product materials pdf (max. 2)to be enslosed in a separate file
- 5. product photo(s) jpg (max. 2) to be enclosed in a separate file

## Clean...

...everywhere in every way

**Covid Products Profile** 



| Product description, area of app  | olication: | Dosage:                  | Temperature/<br>Contact time:  | Packaging:  |
|---|------------|--------------------------|--|---|
| Disinfecting hand cleaning  |            |                          |  |   |
| INNO-SEPT bactericidal (MRSA), fungicidal, tuberculoci<br>(encapsulated viruses incl. H5N1), Hepatitis B virus, disi<br>Cleaning and disinfection in one step.  |            | use undiluted, 3 ml      | B,Y,T: 1 minute<br>s.V: 2 minutes<br>F: 5 minutes<br>at room temperature | 500 ml pump bottle<br>1/5/20 litre canister             |
| LILLA-SEPT bactericidal (MRSA), fungicidal, tuberculocidal, selective virucidal<br>(encapsulated viruses incl. H5N1), Hepatitis B virus, disinfectant hand cleaner.<br>Cleaning and disinfection in one step. |            | use undiluted, 3 ml      | B,Y,T: 1 minute<br>s.V: 2 minutes<br>F: 5 minutes<br>at room temperature | 500 ml pump bottle<br>1/5/20 litre canister             |
| Hand disinfection   |            |                          |  |   |
| INNO-SEPT GÉL EXTRA alcohol based hand disinfects<br>(MRSA), fungicidal, tuberculocidal, virucidal (ind. H1N1)  |            | use undiluted,<br>3-5 ml | 30 seconds/<br>or drying   | 100/500 ml flacon,<br>1litre bottle<br>5 litre canister |
| INNO-SEPT FRESH DISINFECTANT SOLUTION alcoh<br>surface disinfectant, bactericide (MRSA), yeasticidal, tub<br>virucidal (encapsulated viruses).  |            | use undiluted, 3 ml      | 30 seconds<br>T: 60 seconds at<br>room temperature                       | 1 litre pump bottle<br>5/20 litre canister              |
| INNO-SEPT FRESH DISINFECTING WIPES hand and bactericide (MRSA), yeasticidal, tuberculocidal, selective viruses)   |            |                          | 30 seconds<br>T: 60 seconds at<br>room temperature                       | 50 wipes/ case  |

B-bactericidal F-fungicidal Y-yeasticidal V-virucidal T-tuberculocidal